Rapid and Specific Screening Assay for *KRAS* Oncogene Mutation by a Novel Gene Amplification Method

AKIRA KUMASAKA¹, NAOYUKI MATSUMOTO², SHOTARO MUKAE², TAIICHI KITANO², HIROYASU NOGUCHI³, HIDERO OHKI³, KAZUO KOMIYAMA² and TOMOHIRO ANDO¹

¹Department of Oral and Maxillofacial Surgery, Tokyo Women's Medical University, School of Medicine, Shinjyuku-ku, Tokyo, Japan; Departments of ²Pathology, and ³Oral and Maxillofacial Surgery, Nihon University School of Dentistry, Chiyoda-ku, Tokyo, Japan

Abstract. Background: Epidermal growth factor receptor (EGFR) is a target of molecular therapeutics for colorectal cancer. However, mutation of Kirsten rat sarcoma viral oncogene homolog (KRAS) gene at codons 12 and 13 attenuates the therapeutic effect of anti-EGFR therapies. Therefore, the detection of KRAS gene mutation is important for therapeutic decision-making. Materials and Methods: KRAS gene mutation at codons 12 (c.34G>T, c.35G>C, c.35G>A) and 13 (c.38G>A) in six cancer cell lines were investigated. A loop-mediated isothermal amplification-based procedure was developed that employed peptide nucleic acid to suppress amplification of the wild-type allele. Results: This mutation-oriented gene-amplification procedure can amplify the DNA fragment of the KRAS gene with codon 12 and codon 13 mutation within 30 min. Moreover, boiled cells can work as template resources. Conclusion: This newly developed procedure can be useful for patient stratification for anti-EGFR therapies.

Epidermal growth factor (EGFR)-1 induces cell proliferation, inhibition of apoptosis, angiogenesis, and invasion *via* Kirsten rat sarcoma viral oncogene homolog (KRAS) activation in cancer cells (1). This makes EGFR1 a target molecule for cancer therapeutics against several types of carcinoma, including colorectal adenocarcinoma, head and neck squamous cell carcinoma, and non-small cell lung carcinoma (2-6). Recently, antibodies against EGFR, and EGFR tyrosine kinase inhibitors have been used in patients

Correspondence to: Naoyuki Matsumoto, Department of Oral and Maxillofacial Surgery, Nihon University School of Dentistry, 1-8-13, Kanda-Surugadai, Chiyoda-ku, Tokyo 101-8310, Japan. E-mail: matsumoto.naoyuki@nihon-u.ac.jp

Key Words: KRAS mutation, LAMP, molecular targeted therapy, stratification biomarker.

with these types of cancers (2-6). Cetuximab and panitumumab are monoclonal antibodies against the extracellular domain of EGFR1 (7, 8) which inhibit the activation of EGFR; these exhibit antitumor activity when used as monotherapy or in combination with chemotherapy or radiation (9). Clinical efficacy with these molecules has been reported for colonic cancer, head and neck cancer, and non-small cell lung carcinoma (10-12). The EGFR tyrosine kinase inhibitors gefitinib and erlotinib inhibit the autophosphorylation of EGFR and exhibit clinical benefit in the treatment of non-small cell lung carcinoma (13, 14).

The *KRAS* gene mutation, however, attenuates the treatment effect of both monoclonal antibodies against EGFR (15, 16) and EGFR tyrosine kinase inhibitors (17). This mutation is the most common abnormality of the dominant oncogenes in human tumors and is a frequent occurrence in the development and progression of adenocarcinomas of the pancreas (90%), colon (50%), thyroid (50%), bladder (50%), and lung (30%) (18).

Mutations of *KRAS* codons 12 or 13 were detected in 37.7% of colorectal adenoma and adenocarcinoma; of these, 80.8% of the specified mutations occurred in codon 12, with 78.1% of all specified mutations occurring at the second base of either codon (19). Various assays for detecting the *KRAS* mutation include polymerase chain reaction (PCR) restriction fragment length polymorphism (RFLP), PCR-single strand conformation polymorphism (SSCP), PCR sequencing, and mutant–allele-specific amplification (MASA); however, all these methods require DNA purification and involve multiple steps (20-24).

The loop-mediated isothermal amplification (LAMP) assay was developed to amplify copies of a particular DNA sequence (25). The LAMP assay is a rapid and highly specific method compared to PCR. LAMP employs a DNA polymerase with high strand-displacement activity and a set of four basic primers that recognize a total of six distinct sequences on the target DNA (26). Recently, the allele-specific LAMP assay was

0250-7005/2016 \$2.00+.40

Table I. Cell lines employed in this study.

Cell line	Origin	KRAS mutation		
MIA PaCa-2	Pancreatic adenocarcinoma	c.34G>T, p.G12C		
RPMI 8226	Myeloma	c.35G>C, p.G12A		
NCI-H1573	Pulmonary adenocarcinoma	c.35G>C, p.G12A		
SNU-1	Gastric carcinoma	c.35G>A, p.G12D		
Γ84	Colorectal adenocarcinoma	c.38G>A, p.G13D		
DMA-MB-231	Breast adenocarcinoma	c.38G>A, p.G13D		
SK-MEL-28	Melanoma	Wild-type		
HSC-3	Oral squamous cell carcinoma	Wild-type		
NCI-H1975	Non-small cell lung carcinoma	Wild-type		
HT-29	Colorectal adenocarcinoma	Wild-type		
Healthy control	Oral mucosal cells from healthy individuals	Wild-type		

Table II. Primer sequences used in this study.

Primer	Sequence	Target
F3	TGGTGGAGTATTTGATAGTGTA	DNA
В3	ATTGTTGGATCATATTCGTCC	DNA
FIP	AGTCATTTTCAGCAGGCCTTATAAT-TTAACCTTATGTGTGACATGTTC	DNA
BIP	TAAACTTGTGGTAGTTGGAGCTG-TGAATTAGCTGTATCGTCAAGG	DNA
LB	GTGGCGTAGGCAAGAGTG	DNA
LP	CCTACGCCACCAGCTCC	PNA

F3; Forward outer primer, B3; back outer primer, FIP; forward inner primer, BIP; back inner primer, LB; loop-back primer, LP; lock primer.

developed to detect Janus kinase 2 (p.V617F) mutation in chronic myeloproliferative neoplasms (27). The allele-specific LAMP assay employs four basic primers, one self-annealed primer and peptide nucleic acid (PNA), which have been shown to be effective for the enrichment of mutant alleles in PCR screening (28). The PNA could not serve as a primer for polymerization because it is not recognized by the polymerase as DNA and, thus, could not be a substrate for exonuclease activities of Taq polymerase (29).

In this study, we report a novel procedure, mutationoriented LAMP assay, based on the principle of allelespecific LAMP assay. This assay used four basic primers with a back-loop primer (30) and a PNA-based lock primer. The mutation-oriented LAMP assay enabled a rapid and specific detection of KRAS mutation at codons 12 and 13.

Materials and Methods

Ethics statement. This study was approved by the Ethics Committee of the Nihon University School of Dentistry (2009-2). Written informed consent from the donor was obtained for the use of the oral mucosal cells in this research.

Sample preparation. The DNA samples for the LAMP assay were prepared from a pancreatic adenocarcinoma cell line (MIA PaCa-2: c.34G>T, p.G12C), a myeloma cell line (RPMI 8226: c.35G>C,

p.G12A), a pulmonary adenocarcinoma cell line (NCI-H1573: c.35G>C, p.G12A), a gastric carcinoma cell line (SNU-1: c.35G>A, p.G12D), a colorectal adenocarcinoma cell line (T84: c.38G>A, p.G13D), and a breast adenocarcinoma cell line (DMA-MB-231: c.38G>A, p.G13D. These cell lines have point-mutations in KRAS codons 12 or 13. For KRAS wild-type samples, the melanoma cell line (SK-MEL-28), oral squamous cell carcinoma cell line (HSC-3), non-small cell lung carcinoma cell line (NCI-H1975), colorectal adenocarcinoma cell line (HT-29), and oral mucosal cells from healthy individuals were utilized (Table I). All cell lines were obtained from the American Type Culture Collection (Manassas, VA, USA), except for HSC-3 (Riken Bio Resource Center, Ibaraki Japan) and cultured in RPMI-1640 (Gibco, Carlsbad, CA, USA) supplemented with 10% fetal bovine serum (Gibco), and maintained under standard cell culture conditions at 37°C and 5% CO2 in a humid environment. Isolated DNA and boiled cells were used as the template for amplification. Genomic DNA from each cell line was isolated using the DNeasy mini kit (Qiagen, Hilden, Germany) according to the manufacturer's instructions. Both the quality and quantity of DNA were assessed by A260/A280 ratio on a NanoDrop-1000 (Thermo Fisher Scientific Inc., Waltham, MA, USA). Moreover, crude cell lysates were used as LAMP templates. Cells were washed and suspended in phosphate-buffered saline and boiled for 1 min to obtain the crude lysate.

Primer set. A set of five primers were designed for LAMP and mutation-oriented LAMP methods: forward and backward outer primers (F3, B3), forward and backward inner primers (FIP, BIP),




Figure 1. Schematic diagram of loop-mediated isothermal amplification (LAMP) and mutation-oriented LAMP method. The LAMP primer set was designed to bracket the Kirsten rat sarcoma viral oncogene homolog (KRAS) codon 12 and 13 region. A lock primer was designed for the wild-type allelic sequence and tightly bound to the sequence of wild-type allele, inhibiting the LAMP reaction, with resultant amplification of the wild-type allele. B3: Backward outer primer; BIP: backward inner primer; F3: forward outer primer; FIP: forward inner primer.

and back-loop primer (loop B). These primers recognize an area of *KRAS*, including codons 12 and 13.

PNA with exact homology to the wild-type allele of *KRAS*, compatible with codons 12 and 13, was used as the lock primer. In addition to the non-recognition by DNA polymerase making PNA ineligible as a primer for polymerization, the melting temperature of a fully matched PNA–DNA duplex is higher than that of DNA–DNA of the same length, and a single mismatch destabilizes the PNA–DNA hybrids, causing a melting temperature shift of 10-18°C (Figure 1) (31). Thus, the amplification of wild-type allele of *KRAS* gene is inhibited in the mutation-oriented LAMP method. These primer sequences are summarized in Table II.

Mutation-oriented LAMP reaction. The LAMP reaction was performed with an Isothermal Master Mix (Nippon Gene Material Co., Ltd., Tokyo, Japan), which contains DNA polymerase, thermostable pyrophosphatase, reaction buffer, dNTPs, and a double-strand DNA-binding fluorescein dye. The reaction mixture was incubated at 61°C for 35 min. The time to fluorescence positive values (Tp) were monitored on a thermally controlled fluorimeter (Genie® II; Optigene Limited, Horsham, West Sussex, UK) and the Tp values of amplification were automatically recorded. The LAMP reaction was carried out in a 25-µl reaction mixture, containing 40 pmol of each FIP and BIP, 5 pmol of each F3 and B3, specified amount of lock primer, 15 µl of Isothermal Master Mix, distilled water, and target DNA. The lock primer was excluded from the conventional LAMP reaction.

For further confirmation, the amplicons were sequenced with the use of a Big Dye Terminator v3.1 Cycle Sequencing kit (Applied Biosystems, Foster city, CA, USA) and a 3130 XL Genetic analyser (Applied Biosystems) according to the manufacturer's instructions. Sequencing primers were designed to be within the LAMP amplicon.

Determination of optimum temperature. To determine the optimum conditions, the temperature and concentration of PNA were examined. The reaction mix was incubated at 60-65°C for 30 min and Tp was monitored. Alternatively, PNA was added to the reaction mix at different concentrations of 5, 10, 15, 16.75, 20, and 40 pmol, respectively.

Statistical analysis. To examine the difference of Tp values, the Kruskal–Wallis H-test was performed. All experiments were performed five times. Differences were considered statistically significant at *p*-values of less than 0.05.

Results

LAMP and mutation-oriented LAMP reaction. We compared the conventional LAMP and mutation-oriented LAMP reaction for the *KRAS* gene. In the LAMP reaction, all examined samples exhibited a positive reaction for *KRAS*. The amplification of the gene was observed in 10-15 min. There was no difference in the Tp value of the *KRAS* wild-type and mutation groups. In contrast, the mutation-oriented LAMP reaction delayed gene amplification in the *KRAS* wild-type group (NCI-H1975, HSC-3, HT-29, and healthy control). In the LAMP reaction, the Tp value was not significantly different for the *KRAS* mutation (15'56) and wild-type groups (14'56). However, in the mutation-oriented LAMP reaction, the Tp value significantly differed for the *KRAS* mutation (18'34) and wild-type groups (35'42) (*p*<0.0001; Figure 2).




Figure 2. Loop-mediated isothermal amplification (LAMP) and the mutation-oriented LAMP reaction in cancer cell lines. The mutation-oriented LAMP reaction was suppressed in the HSC-3 cell line (oral squamous cell carcinoma), NCI-H1975 cell line (lung adenocarcinoma), and healthy oral mucosal cells, which have the wild-type allele for Kirsten rat sarcoma viral oncogene homolog (KRAS) (solid line). The LAMP reaction amplifies all KRAS alleles (broken line).

We measured the difference between the Tp value of mutation-oriented LAMP reaction and the LAMP reaction. The differences were: 1'00 in T84, -0'12 in NCI-H1573, 2'54 in MIA PaCa-2, 1'00 in SNU-1, 14'51 in HSC-3, 10'14 in NCI-H1975, 12'08 in HT-29, and 13'33 in healthy controls (Table III).

The DNA sequences of each amplicon were evaluated and were found to be compatible with the estimated sequences. *Optimization of reaction conditions*. To explore optimal conditions, the reaction temperature and lock primer concentration were examined. The time at which the upstroke of fluorescence level elevation was recorded, and the




Figure 3. Detection limit for Kirsten rat sarcoma viral oncogene homolog (KRAS) mutation by loop-mediated isothermal amplification (LAMP) reaction (A) and mutation-oriented LAMP reaction (B). Both LAMP and mutation-oriented LAMP reactions amplify the KRAS gene at concentrations over 10 pg; however, the time to positive fluorescence (Tp) value was delayed at concentrations of 10 and 100 pg in the mutation-oriented LAMP reaction.

Table III. Differences between loop-mediated isothermal amplification (LAMP) and mutation-oriented LAMP reaction.

	Mutant allele				Wild allele			
	T84	NCI-H 1573	MIA PaCa-2	SNU-1	NCI-H 1975	HSC-3	HT-29	Healthy control
LAMP	13'57	16'46	16'08	17'43	14'56	15'52	14'38	13'11
Mutation-oriented LAMP	14'57	15'58	19'02	18'43	35'10	30'43	26'46	26'44
$\Delta \mathrm{Tp}$	1'00	-0'12	2'54	1'00	20'14	14'51	12'08	13'33
p-Value	>0.05	< 0.05	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01

ΔTp: Time to positive value (Tp) of mutation-oriented LAMP minus the Tp value of LAMP.

difference of Tp value between the KRAS wild-type and mutation groups was examined. To compare the difference of reaction temperature, reaction mixtures were incubated at 60-65°C. Differences of Tp value in the two groups at different reaction temperatures are presented in Table IV. At 65°C, the Tp value was faster and the difference between the wild and mutant types was the greatest (Table IV). To compare the difference of the concentrations of lock primer, we conducted a mutation-oriented LAMP reaction with different PNA concentrations. Differences of Tp value between the *KRAS* wild-type (healthy control) and *KRAS* mutation type (T84) are presented in Table V. There was little difference beyond the concentration of 16.75 pmol. Therefore, we conducted both LAMP and mutation-oriented LAMP reaction at 61°C with 16.75 pmol PNA for further experiments (Table V).

To examine the detection limit of the mutation-oriented LAMP reaction, a template gDNA of T-84 cells was added into the reaction mix at concentrations of 100 fg to 15 ng.

Both mutation-oriented LAMP and conventional LAMP reactions were conducted, and the results showed that both reaction procedures amplified the mutation allele of the *KRAS* gene over 10 pg of template gDNA; however, the Tp value was delayed at concentrations of 10 and 100 pg in the mutation-oriented LAMP reaction (Figure 3). Therefore, we infer that the mutation-oriented LAMP reaction requires a template gDNA of mutation allele of over 1 ng.

Furthermore, we performed the assay using genomic DNA from T-84 cell line gradually diluted (100, 10, 1, 0.1, 0.01, 0.001 and 0%) into wild-type DNA from healthy controls. The reproducible sensitivity of mutation-oriented LAMP was 10%, detected in all five repeated examinations (100%). The 1% dilution was detected in 60% (Table VI).

Mutation-oriented LAMP analysis directly from crude samples. We also investigated boiled cell lines as template. Each boiled cell works as the template for LAMP and




Figure 4. Detection of Kirsten rat sarcoma viral oncogene homolog (KRAS) mutation using crude samples. Crude samples were obtained by washing and boiling several cultured cells in phosphate-buffered saline. Mutated KRAS codons 12 and 13 were detected from over 100 cells in the reaction mixture (500 cells; solid line, 100 cells; gray line). Moreover, mutant alleles could be detected from over 10 cells (broken line) in SNU-1 cell.

mutation-oriented LAMP systems. We examined the detection limit of number of cells, and amplification was observed in 500 cells (T84: 16'34±0'11, NCI-H1573: 18'03±0'07, MIA PaCa-2: 17'48±0'19, SNU-1: 15'08±0'14) and 100 cells (T84: 25'07±0'16, NCI-H1573: 23'07±0'09, MIA PaCa-2: 22'39±0'18, SNU-1: 15'51±0'23); however, amplification could not be detected in samples of 10 cells, except for SNU-1 cells (16'48) (Figure 4).

Discussion

The detection of the *KRAS* mutation at codons 12 and 13 is quite important for cancer therapeutics using antibodies against EGFR and EGFR tyrosine kinase inhibitor. We successfully developed a mutation-oriented LAMP assay using several cultured cell lines. This assay was able to amplify the DNA fragment of the mutated *KRAS* gene within 20 min. Lock primers made of PNA can inhibit gene amplification in *KRAS* wild-type codons 12 and 13. Thus, the assay is able to amplify all possible mutations in *KRAS* codons 12 and 13.

Various assays for KRAS mutation detection exist, including PCR-RFLP, PCR-SSCP, MASA, and other PCRbased technologies or DNA sequencing (20-24). PCR-RFLP requires DNA amplification, restriction enzyme digestion, and electrophoresis. PCR-SSCP requires DNA amplification, denaturation of amplicons, and polyacrylamide gel electrophoresis or capillary-gel electrophoresis. MASA requires DNA amplification with several PCR primer sets for each mutation as well as electrophoresis. All of these PCRbased methods are time consuming (shortest time, 3.5 h), require careful DNA extraction, and involve several steps (DNA extraction, PCR, electrophoresis, etc.) (32). DNA sequencing can detect unknown mutations; however, this is a time- and cost-ineffective method. In contrast, the mutationoriented LAMP reaction enables identification of gene mutation within 20 minutes.

The mutation-oriented LAMP method enables targeting of all *KRAS* mutations in codons 12 and 13 at once. Jakovljevic *et al.* reported that 34.7% of patients with metastatic colorectal cancer have *KRAS* point mutations in codons 12

Table IV. Differences in mutation-oriented loop-mediated isothermal amplification (LAMP) reaction speed based on temperature.

	60°C	61°C	62°C	63°C	64°C	65°C
Mutant allele						
T84	18'04±0'42	14'57±0'18	14'41±0'21	13'32±0'19	13'25±0'15	12'55±0'17
NCI-H1573	18'49±0'11	19'12±1'02	16'41±0'16	14'17±0'12	13'25±0'24	13'40±0'14
Mutant allele average	18'27±0'51	17'34±1'11	15'49±0'31	13'54±0'24	13'25±0'33	13'16±0'51
Normal allele						
NCI-H1975	37'34±2'54	35'42±3'07	29'11±1'11	28'02±1'06	22'10±0'58	19'25±0'51
Healthy controls	28'04±1'01	35'57±3'01	22'26±1'07	19'47±0'37	18'55±0'12	18'55±1'03
Normal allele average	32'48±3'24	35'19±3'11	25'57±2'13	23'54±1'56	20'33±0'58	19'25±0'51
Tp value of mutant average – normal allele average	14'21±3'51	18'45±3'34	10'08±3'11	10'00±2'35	7'08±2'47	6'53±3'15
<i>p</i> -Value	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01

Tp: Time to positive fluorescence.

Table V. Differences in mutation-oriented loop-mediated isothermal amplification (LAMP) reaction based on lock primer concentration.

	Lock primer concentration (pmol/reaction)						
	5	10	15	16.75	20	40	
T84	14'20±0'11	14'55±0'13	15'05±0'19	16'54±0'23	16'10±0'21	17'40±0'33	
Healthy controls	15'50±0'21	18'25±0'32	22'20±1'47	23'50±2'05	22'10±1'43	26'25±2'01	
Tp for T84 –Tp for healthy control <i>p</i> -Value	01'31±0'41 >0.05	03'30±2'43 <0.05	07'15±4'21 <0.01	06'46±5'10 <0.01	06'00±4'51 <0.01	08'45±5'23 <0.01	

Tp: Time to positive fluorescence.

Table VI. Sensitivity of mutation-oriented loop-mediated isothermal amplification (LAMP) assay in mixed-cell populations.

	Mutant allele concentration (%)						
	100	10	1	10^{-1}	10-2	10-3	0
LAMP	14'24±0'07	16'27±0'19	22'54±0'33	24'19±0'47	26'33±0'51	28'45±1'09	28'42±1'11
Mutation-oriented LAMP	14'40±0'05	14'48±0'11	15'08±0'13	15'34±0'17	15'12±0'21	15'06±0'20	15'08±0'22
ΔTp	-0'04±0'06	1'39±0'17	7'42±0'29	8'48±0'33	11'21±0'45	13'38±0'46	13'34±0'47
<i>p</i> -Value			< 0.01	< 0.01	< 0.01	< 0.01	< 0.01

ΔTp:Time to positive fluorescence (Tp) value of mutation-oriented LAMP minus Tp value of LAMP.

and 13. Among them, the frequencies of mutations in codon 12 were: p.G12D, 43.9%; p.G12V, 21.2%; p.G12A, 10.6%; p.G12C, 7.6%; p.G12S, 4.5%; p.G12R, 1.5%; the frequency of the only mutation in codon 13 was 10.6% for p.G13D (33). In the mutation-oriented LAMP method, the wild-type allele of *KRAS* was masked by lock primer, composed of PNA with a complementary sequence for *KRAS* codons 12 and 13. PNAs are DNA mimics, wherein the phosphoribose backbone is replaced by a peptide-like repeat of (2-aminoethyl)-glycine

units (34). Due to this chemical difference, PNA differs from DNA molecules in several aspects. PNA–DNA hybrids have a higher thermal stability compared to the corresponding DNA–DNA hybrids. PNA–DNA hybrids are more easily destabilized by single base-pair mismatches than the corresponding DNA–DNA hybrids. Thus, PNAs cannot serve as primer molecules in PCR (31, 35). Therefore, if there are mutations in codons 12 and 13, the lock primer detaches from the template DNA, and the LAMP reaction is initiated.

Another advantage of this novel LAMP method was its detection ability in the crude sample. Conventional *KRAS* mutation detection methods require DNA extraction from samples. However, we were able to detect the *KRAS* mutations without performing DNA extraction. This was due to the high sensitivity of the LAMP method (26). In this study, we successfully detected the *KRAS* mutation from 100 boiled cells in T84 (c.38G>A), NCI-H-1573 (c.35G>C), and MIA PaCa-2 (c.34G>T) cell lines. Moreover, we detected mutations from as few as 10 cells in the SNU-1 (c.35G>A) cell line. These differences could be attributable to the differences of the substituted base and its position.

In summary, the mutation-oriented LAMP method detects *KRAS* gene mutation in codons 12 and 13 by using lock primers. As this analysis can be carried out without DNA extraction, this method could be applied for clinical examination in future applications.

Acknowledgements

This study was supported by Grant-in-Aid for Young Scientists (B) to Shotaro Mukae (No.24792011), and Grant-in-Aid for Scientific Research (C) to Naoyuki Matsumoto from the Ministry of Education, Cultures, Sports, Science and Technology, Japan.

References

- 1 Raymond E, Faivre S and Armand JP: Epidermal growth factor receptor tyrosine kinase as a target for anticancer therapy. Drugs 60(Suppl 1): 15-23, 2000.
- 2 Di Gennaro E, Barbarino M, Bruzzese F, De Lorenzo S, Caraglia M, Abbruzzese A, Avallone A, Comella P, Caponigro F, Pepe S and Budillon: A Critical role of both p27^{KIP1} and p21^{CIP1/WAF1} in the antiproliferative effect of ZD1839 ('Iressa'), an epidermal growth factor receptor tyrosine kinase inhibitor, in head and neck squamous carcinoma cells. J Cell Physiol 195: 139-150, 2003.
- 3 Cunningham D, Humblet Y, Siena S, Khayat D, Bleiberg H, Santoro A, Bets D, Mueser M, Harstrick A, Verslype C, Chau I and Van Cutsem E: Cetuximab monotherapy and cetuximab plus irinotecan in irinotecan-refractory metastatic colorectal cancer. N Engl J Med 351: 337-345, 2004.
- 4 Burtness B, Goldwasser MA, Flood W, Mattar B and Forastiere AA; Eastern Cooperative Oncology Group: Phase III randomized trial of cisplatin plus placebo compared with cisplatin plus cetuximab in metastatic/recurrent head and neck cancer: an Eastern Cooperative Oncology Group study. J Clin Oncol 23: 8646-8654, 2005.
- 5 Linardou H, Dahabreh IJ, Kanaloupiti D, Siannis F, Bafaloukos D, Kosmidis P, Papadimitriou CA and Murray S: Assessment of somatic K-RAS mutations as a mechanism associated with resistance to EGFR-targeted agents: a systematic review and meta-analysis of studies in advanced non-small-cell lung cancer and metastatic colorectal cancer. Lancet Oncol 9: 962-972, 2008.
- 6 Pirker R, Pereira JR, Szczesna A, von Pawel J, Krzakowski M, Ramlau R, Vynnychenko I, Park K, Yu CT, Ganul V, Roh JK, Bajetta E, O'Byrne K, de Marinis F, Eberhardt W, Goddemeier T, Emig M and Gatzemeier U; FLEX Study Team: Cetuximab plus

- chemotherapy in patients with advanced non-small-cell lung cancer (FLEX): an open-label randomised phase III trial. Lancet *373*: 1525-1531, 2009.
- 7 Mendelsohn J and Baselga: J The EGF receptor family as targets for cancer therapy. Oncogene *19*: 6550-6565, 2000.
- 8 Yang XD, Jia XC, Corvalan JR, Wang P and Davis CG: Development of ABX-EGF, a fully human anti-EGF receptor monoclonal antibody, for cancer therapy. Crit Rev Oncol Hematol 38: 17-23, 2001.
- 9 Herreros-Villanueva M, Chen CC, Yuan SS, Liu TC and Er TK: KRAS mutations: analytical considerations. Clin Chim Acta 431: 211-220, 2014.
- 10 Hanna N, Lilenbaum R, Ansari R, Lynch T, Govindan R, Jänne PA and Bonomi P: Phase II trial of cetuximab in patients with previously treated non-small-cell lung cancer. J Clin Oncol 24: 5253-5258, 2006.
- 11 Bonner JA, Harari PM, Giralt J, Azarnia N, Shin DM, Cohen RB, Jones CU, Sur R, Raben D, Jassem J, Ove R, Kies MS, Baselga J, Youssoufian H, Amellal N, Rowinsky EK and Ang KK: Radiotherapy plus cetuximab for squamous-cell carcinoma of the head and neck. N Engl J Med 354: 567-578, 2006.
- 12 Saltz LB, Meropol NJ, Loehrer PJ Sr, Needle MN, Kopit J and Mayer RJ: Phase II trial of cetuximab in patients with refractory colorectal cancer that expresses the epidermal growth factor receptor. J Clin Oncol 22: 1201-1208, 2004.
- 13 Pérez-Soler R, Chachoua A, Hammond LA, Rowinsky EK, Huberman M, Karp D, Rigas J, Clark GM, Santabárbara P and Bonomi P: Determinants of tumor response and survival with erlotinib in patients with non-small-cell lung cancer. J Clin Oncol 22: 3238-3247, 2004.
- 14 Kim KS, Jeong JY, Kim YC, Na KJ, Kim YH, Ahn SJ, Baek SM, Park CS, Park CM, Kim YI, Lim SC and Park KO: Predictors of the response to gefitinib in refractory non-small cell lung cancer. Clin Cancer Res 11: 2244-2251, 2005.
- 15 Lièvre A, Bachet JB, Le Corre D, Boige V, Landi B, Emile JF, Côté JF, Tomasic G, Penna C, Ducreux M, Rougier P, Penault-Llorca F and Laurent-Puig P: KRAS mutation status is predictive of response to cetuximab therapy in colorectal cancer. Cancer Res 66: 3992-3995, 2006.
- 16 Amado RG, Wolf M, Peeters M, Van Cutsem E, Siena S, Freeman DJ, Juan T, Sikorski R, Suggs S, Radinsky R, Patterson SD and Chang DD: Wild-type *KRAS* is required for panitumumab efficacy in patients with metastatic colorectal cancer. J Clin Oncol 26: 1626-1634, 2008.
- 17 Pao W, Wang TY, Riely GJ, Miller VA, Pan Q, Ladanyi M, Zakowski MF, Heelan RT, Kris MG and Varmus HE: *KRAS* mutations and primary resistance of lung adenocarcinomas to gefitinib or erlotinib. PLoS Med 2: e17, 2005.
- 18 Bissada E, Abboud O, Abou Chacra Z, Guertin L, Weng X, Nguyen-Tan PF, Tabet JC, Thibaudeau E, Lambert L, Audet ML, Fortin B and Soulières D: Prevalence of *K-RAS* codons 12 and 13 mutations in locally advanced head and neck squamous cell carcinoma and impact on clinical outcomes. Int J Otolaryngol 2013: 848021, 2013.
- 19 Andreyev HJ, Norman AR, Cunningham D, Oates JR and Clarke PA: Kirsten ras mutations in patients with colorectal cancer: the multicenter "RASCAL" study. J Natl Cancer Inst 90: 675-684, 1998.
- 20 Hatzaki A, Razi E, Anagnostopoulou K, Iliadis K, Kodaxis A, Papaioannou D, Labropoulos S, Vasilaki M, Kosmidis P, Saetta

- A, Mihalatos M and Nasioulas G: A modified mutagenic PCR-RFLP method for *K-RAS* codon 12 and 13 mutations detection in NSCLC patients. Mol Cell Probes *15*: 243-247, 2001.
- 21 Mulcahy HE, Lyautey J, Lederrey C, qi Chen X, Anker P, Alstead EM, Ballinger A, Farthing MJ and Stroun M: A prospective study of *K-RAS* mutations in the plasma of pancreatic cancer patients. Clin Cancer Res 4: 271-275, 1998.
- 22 Yamashita K, Tatebayashi T, Shinoda H and Okayasu I: Simplified rapid non-radioactive PCR-SSCP method applied to K-RAS mutation analysis. Pathol Int 46: 801-804, 1996.
- 23 Takeda S, Ichii S and Nakamura Y: Detection of K-RAS mutation in sputum by mutant-allele-specific amplification (MASA). Hum Mutat 2: 112-117, 1993.
- 24 Mok SC, Bell DA, Knapp RC, Fishbaugh PM, Welch WR, Muto MG, Berkowitz RS and Tsao SW: Mutation of *K-RAS* protooncogene in human ovarian epithelial tumors of borderline malignancy. Cancer Res 53: 1489-1492, 1993.
- 25 Iwamoto T, Sonobe T and Hayashi K: Loop-mediated isothermal amplification for direct detection of Mycobacterium tuberculosis complex, M. avium, and M. intracellulare in sputum samples. J Clin Microbiol 41: 2616-2622, 2003.
- 26 Notomi T, Okayama H, Masubuchi H, Yonekawa T, Watanabe K, Amino N and Hase T: Loop-mediated isothermal amplification of DNA. Nucleic Acids Res 28: E63, 2000.
- 27 Minnucci G, Amicarelli G, Salmoiraghi S, Spinelli O, Guinea Montalvo ML, Giussani U, Adlerstein D and Rambaldi A: A novel, highly sensitive and rapid allele-specific loop-mediated amplification assay for the detection of the *JAK2V617F* mutation in chronic myeloproliferative neoplasms. Haematologica 97: 1394-1400, 2012.
- 28 Demers DB, Curry ET, Egholm M and Sozer AC: Enhanced PCR amplification of *VNTR* locus D1S80 using peptide nucleic acid (PNA). Nucleic Acids Res 23: 3050-3055, 1995.
- 29 Tatsumi K, Mitani Y, Watanabe J, Takakura H, Hoshi K, Kawai Y, Kikuchi T, Kogo Y, Oguchi-Katayama A, Tomaru Y, Kanamori H, Baba M, Ishidao T, Usui K, Itoh M, Cizdziel PE, Lezhava A, Ueda M, Ichikawa Y, Endo I, Togo S, Shimada H and Hayashizaki Y: Rapid screening assay for KRAS mutations by the modified smart amplification process. J Mol Diagn 10: 520-526, 2008.

- 30 Nagamine K, Hase T and Notomi T: Accelerated reaction by loop-mediated isothermal amplification using loop primers. Mol Cell Probes 16: 223-229, 2002.
- 31 Kyger EM, Krevolin MD and Powell MJ: Detection of the hereditary hemochromatosis gene mutation by real-time fluorescence polymerase chain reaction and peptide nucleic acid clamping. Anal Biochem *260*: 142-148, 1998.
- 32 Hoshi K, Takakura H, Mitani Y, Tatsumi K, Momiyama N, Ichikawa Y, Togo S, Miyagi T, Kawai Y, Kogo Y, Kikuchi T, Kato C, Arakawa T, Uno S, Cizdziel PE, Lezhava A, Ogawa N, Hayashizaki Y and Shimada H: Rapid detection of epidermal growth factor receptor mutations in lung cancer by the SMart-Amplification Process. Clin Cancer Res 13: 4974-4983, 2007.
- 33 Jakovljevic K, Malisic E, Cavic M, Krivokuca A, Dobricic J and Jankovic R: KRAS and BRAF mutations in Serbian patients with colorectal cancer. J BUON 17: 575-580, 2012.
- 34 Nielsen PE, Egholm M, Berg RH and Buchardt O: Sequence-selective recognition of DNA by strand displacement with a thymine-substituted polyamide. Science 254: 1497-1500, 1991.
- 35 Thiede C, Bayerdörffer E, Blasczyk R, Wittig B and Neubauer A: Simple and sensitive detection of mutations in the RAS protooncogenes using PNA-mediated PCR clamping. Nucleic Acids Res 24: 983-984, 1996.

Received January 20, 2016 Revised March 8, 2016 Accepted March 17, 2016